



# HemCon<sup>®</sup> OneStop<sup>™</sup> Vascular

HemCon<sup>®</sup> OneStop<sup>™</sup> Vascular is a hemostatic dressing for the external, temporary control of severely bleeding wounds intended for emergency use.


HemCon<sup>®</sup> OneStop<sup>™</sup> Vascular also controls bleeding in patients following hemodialysis and is indicated for the control of bleeding from the skin at percutaneous needle access, vascular access and percutaneous catheter access sites.

HemCon<sup>®</sup> OneStop<sup>™</sup> Vascular is fabricated from chitosan, a naturally occurring, biocompatible polysaccharide.

## Features:

An ideal solution for safe bleeding control.


- Rapid control of bleeding
- Hemostasis in patients on anticoagulation therapy
- Provides an antibacterial barrier to MRSA, VRE, and many others<sup>1</sup>
- Maintains structural Integrity - won't crack, crumble, shed, or become saturated
- Clinically proven to reduce hold times
- Can be used as an adjunct with mechanical closure devices



<sup>1</sup> Data on file at Tricol, in vitro study.  
Testing performed by an independent, certified, contract laboratory

## Easy Application

- 1 Allow nickel-sized (~2 cm) drop of blood to form at puncture site. Blood is required for dressing to adhere here. (Do not cleanse puncture site or moisten with saline solution)
- 2 With printed side facing up, place dressing directly on puncture site. Dressing can be cut to size. Do not remove the backing
- 3 Hold digital pressure until bleeding is controlled
- 4 After bleeding has stopped, secure OneStop™ Vascular with appropriate dressing (not included)



## Easy Removal

- 1 Recheck the wound for potential bleeding as necessary. If hemostasis is not achieved or for recurrent bleeding, remove dressing with saline or water and re-apply a new dressing until hemostasis is achieved
- 2 Remove dressing within 48 hours by irrigating with saline or water while gently pulling up on the corner of the dressing

## Intended Use / Indications:

The dressing is intended for the local management of bleeding wounds and to provide a barrier to bacterial penetration of the dressing in all patients and for the promotion of rapid control (hemostasis) of bleeding in patients following hemodialysis and for those on anticoagulation therapy.

The dressing is indicated for the following wounds: skin surface puncture sites, vascular procedure sites, and sites involving percutaneous catheters.

## Unique Chitosan-based Biotechnology

OneStop™ products are composed of chitosan, a naturally occurring, biocompatible polysaccharide derived from shrimp shells.

Chitosan's positive molecular charge attracts negatively charged red blood cells, similar to a magnet. As the red blood cells are drawn to the bandage, a clot is formed in the wound.

The result:

- A tight seal over the wound site
- Fast hemostasis - separate from, and supportive of, the body's natural ability to clot
- An antibacterial barrier<sup>1</sup>

## Ordering Information

Ref no.	Description	Qty/box
1102	HemCon® OneStop™ 4 cm x 4 cm	10
1103	HemCon® OneStop™ 5 cm x 5 cm	10

### Contact Details, Nordics

customer@mermaidmedical.com  
mermaidmedical.com

### Contact Details, NL

info@mermaidmedical.nl  
mermaidmedical.com