




FlexGRIP® Catheter Securement Device



The FlexGRIP® device is suitable for securing and stabilizing catheters from the M-Drain® series, as well as other catheters ranging from 5-18 French. It is made of soft materials and provides an easy-to-use solution.

The soft, low profile retainer and pre-mounted fixation leash securely holds the catheter in place, with an additional strap that locks the soft clip to avoid accidental loosening. The adhesive is skin-safe and hypoallergenic, designed to prevent MARS (medical adhesive-related skin injury) for optimal patient comfort.

FlexGRIP® comes in 2 sizes to ensure a perfect fit for every catheter: 5-11Fr and 11 to 18Fr.

Features:

- Air pillow technology ensures tight catheter fixation without impacting flow rate
- Hypoallergenic and skin-safe adhesive designed to prevent MARS
- Pre-mounted fixation leash for ease of use
- Soft materials and low profile for enhanced patient comfort

Ordering Information

Ref no.	Description	Qty/box
UN511	Securement Device For Percutaneous Catheters · 5Fr-11Fr	50
UN1118	Securement Device For Percutaneous Catheters · 11Fr-18Fr	25

Contact Details

Frydensbergvej 25 · DK-3660 Stenløse
customer@mermaidmedical.com · mermaidmedical.com