



ENDOROTOR



The EndoRotor® is a next generation interventional resection device designed to address the challenges of pancreatic necrosectomy procedures, sessile lesions and incomplete resections.

The EndoRotor® is indicated for endoscopic removal of necrotic pancreas tissue and resection of tissue from within the GI tract.

The tool resects and aspirates tissue simultaneously through a catheter inserted into the working channel of an endoscope.

This interventional tool replaces multiple instruments and the time necessary for instrument exchange.

Features:

- This interventional tool replaces multiple instruments and reduces the time necessary for instrument exchange.
- The necrosectomy tool can overcome the challenge of removing the varied types of necrotic tissue with just one device and without leaving the pancreatic cavity.
- The mucosectomy tool is designed to address the challenge of tissue bridges as well as tissue around the lateral margin. It is designed to address sessile/flat lesions by suctioning and cutting in one motion
- Preloaded applicator
- Damage protection



Necrosectomy catheter

Simultaneously resects and removes necrotic tissue and fluid from a pancreatic cyst through a catheter inserted into the working channel of an endoscope.



Mucosectomy catheter

Simultaneously resects and removes mucosal tissue through a catheter inserted into the working channel of an endoscope.

Ordering Information

EndoRotor® Resection Device

Ref. no.	Description	Qty/box
3.2-EPR-COL-OP	EndoRotor® Catheter Colon mucosectomy for Olympus & Pentax	1
3.2-EPR-EGD-F	EndoRotor® Catheter Gastro mucosectomy for FujiFilm	1
3.2-EPR-EGD-OP	EndoRotor® Catheter Gastro mucosectomy for Olympus & Pentax	1
3.2-PED-EGD-F	EndoRotor® Catheter Necrosectomy for FujiFilm	1
3.2-PED-EGD-OP	EndoRotor® Catheter Necrosectomy for Olympus & Pentax	1
6.0-PED-EGD	EndoRotor®-6.0mm Catheter Necrosectomy for Olympus	1
CATHGUIDE	EndoRotor® Catheter Guide for 6.0-PED-EGD	1
PED-PURGE	Necrosectomy purge kit	10
EPR-TRAP	EndoRotor® Specimen Trap	10
ESS- CONSOLE	EndoRotor® Console	1
ESS-KIT	EndoRotor® Console, roll stand and vacuum pump	1