



# MedOne Ultra

## Programmable automatic biopsy system



MedOne Ultra programmable automatic biopsy system is designed with a powerful spring firing mechanism to obtain a high-quality biopsy sample.

Its one-handed usage and lightweight design make it easy to manipulate. It holds a 25 mm penetration depth and a 20 mm sample notch.


MedOne Ultra features automatic or sequential firing modes;

**D:** Delay firing mode to visualize needle stylet sample notch placement within a lesion before activating the cutting cannula.

**A:** Automatic firing mode releases the stylet and cannula in rapid sequence capturing specimen by pushing a single button.

Features:

- Automatic or sequential firing modes
- 25 mm penetration depth to gain high-quality 20 mm tissue samples
- One-hand use leaves the physician able to operate the ultrasound probe



## Ordering Information

### Automatic Disposable Biopsy System

Ref no.	Gauge	Length	Qty/box
MU14100-00	14	10 CM	5
MU14160-00	14	16 CM	5
MU16100-00	16	10 CM	5
MU16160-00	16	16 CM	5
MU18100-00	18	10 CM	5
MU18160-00	18	16 CM	5
MU18200-00	18	20 CM	5
MU18250-00	18	25 CM	5
MU20100-00	20	10 CM	5
MU20160-00	20	16 CM	5
MU20200-00	20	20 CM	5

### Coaxial Needles Sold Separately

Ref no.	Gauge	Length	Qty/box
CT14100-00	13	5.5 CM	20
CT14160-00	13	11.5 CM	20
CT16100-00	15	5.8 CM	20
CT16160-00	15	11.8 CM	20
CT18100-00	17	5.9 CM	20
CT18160-00	17	11.9 CM	20
CT18200-00	17	15.9 CM	20
CT18250-00	17	20.9 CM	20
CT20100-00	19	5.75 CM	20
CT20160-00	19	11.75 CM	20
CT20200-00	19	15.75 CM	20

### Automatic Disposable Biopsy System with Coaxial Needle

Ref no.	Gauge	Length	Qty/box
MU14100-C0	14	10 CM	5
MU14160-C0	14	16 CM	5
MU16100-C0	16	10 CM	5
MU16160-C0	16	16 CM	5
MU18100-C0	18	10 CM	5
MU18160-C0	18	16 CM	5
MU18200-C0	18	20 CM	5
MU18250-C0	18	25 CM	5
MU20100-C0	20	10 CM	5
MU20160-C0	20	16 CM	5
MU20200-C0	20	20 CM	5

#### Contact Details, Nordics

Frydensbergvej 25 · DK-3660 Stenløse  
[customer@mermaidmedical.com](mailto:customer@mermaidmedical.com) · [mermaidmedical.com](http://mermaidmedical.com)

#### Contact Details, UK

3A Cable Street, Formby · Merseyside L37 3LU  
[customer@mermaidmedical.co.uk](mailto:customer@mermaidmedical.co.uk) · [mermaidmedical.com](http://mermaidmedical.com)