



## Blue Eye™ Submucosal Injection Agent



Blue Eye™ is a medical device that is an injectable liquid composition intended for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early-stage cancers, or other gastrointestinal mucosal lesions, prior to removal with a snare or endoscopic device.

It facilitates endoscopic resection procedures during endoscopic examinations in the upper and lower gastrointestinal tract, such as the esophagus, the stomach.

### Features:

- Ready to use and easy to connect to the injection needle, which shortens the procedure time and minimizes the risk of cross-contamination.
- Has a mucosal elevation time twice as long as a saline solution, which shortens the procedure time, and this reduces the risk of complications.
- Contains methylene blue for visualizing the lesion.
- Contains hyaluronic acid, a biopolymer that occurs naturally in human skin and connective tissue, making Blue Eye™ biocompatible.

### Ordering Information

Ref. no.	Description	Qty/box
TS-905	Blue Eye™ Submucosal Injection Agent (5ml)	10



info@mermaidmedical.nl  
mermaidmedical.com