

# iVascular<sup>®</sup>

therapies for living

## sergeant

peripheral support catheter



To lead complex lesions procedures



[www.ivascular.global](http://www.ivascular.global)

# sergeant

Designed to ease the complex and distal lesions procedures

## Support to cross all lesions

Developed to ensure an outstanding guidewire support:

- Braided catheter with exceptional pushability and torque capacity
- Tapered shaft design to increase support



## Able to go through tortuous anatomy

HYDRAX +, proprietary hydrophilic coating that increases the catheter trackability

## Perfect visibility


The 4 radiopaque markers assess accurately:

- Catheter position
- Lesion length

# sergeant features

- > **Over the wire catheter (OTW)**
- > Tip: **Straight or 30° Angled**


- > **4 tungsten radiopaque markers**
- > **4F** introducer compatible



Product with CE mark, certified by Notified Body 0123

Guidewire compatibility	Effective catheter length (cm)				Tip shape
	65	90	130	150	
0.014"	DSP C14 065 001	DSP C14 090 001	DSP C14 130 001	DSP C14 150 001	Straight
	DSP C14 065 002	DSP C14 090 002	DSP C14 130 002	DSP C14 150 002	Angled
0.018"	DSP C18 065 001	DSP C18 090 001	DSP C18 130 001	DSP C18 150 001	Straight
	DSP C18 065 002	DSP C18 090 002	DSP C18 130 002	DSP C18 150 002	Angled
0.035"	DSP C35 065 001	DSP C35 090 001	DSP C35 130 001	DSP C35 150 001	Straight
	DSP C35 065 002	DSP C35 090 002	DSP C35 130 002	DSP C35 150 002	Angled

Distributed by:



Manufactured by:

Life Vascular Devices Biotech, S.L.  
[www.ivascular.global](http://www.ivascular.global)  
[info@ivascular.global](mailto:info@ivascular.global)