



mermaid medical®

Angel® Catheter

Inferior vena cava filter and
central venous catheter system




 ANGEL®
CATHETER



Angel® Catheter IVC Filter Retrieval, Finally Guaranteed

The Angel® Catheter is an IVC filter permanently attached to a triple lumen central venous catheter.

The smart, elegant design of the Angel® Catheter provides critical care physicians and their patients with prophylactic PE protection and guarantees IVC filter retrieval, a game changing promise that no other IVC filter can make.



As the first and only IVC filter to receive FDA clearance for a prophylactic use indication, the Angel® Catheter provides immediate and effective PE protection during a critical and frequently untreated window for critically ill patients.

During this unprotected time when patients are considered high risk for VTE and contraindicated to anticoagulants, the Angel® Catheter allows for bedside placement and provides newfound peace of mind for clinicians and their patients.

Indications For Angel® Catheter Placement

- The Angel® Catheter is intended for short-term use for the prevention of clinically significant pulmonary embolism (PE) in critically ill patients at high risk for PE or recurrent PE, and recognized contraindications to standard pharmacological thromboprophylaxis therapy.
- The Angel® Catheter is also intended to provide access to the central venous system.

Technical information

Guidewire compatibility	0.035"
Filter size	30 x 50 mm
Catheter profile and useable length	9F, 30cm
Power injection	Distal and Proximal port
MR conditional	1.5 and 3 Tesla





mermaid medical[®]

Mermaid Medical, Inc.

7340 South Alton Way, Suite 11-H
Centennial, CO, 80112, USA
Customercare@mermaidmedical.com
Phone: 720-420-9999

Angel[®] Catheter

Item nr.	Description	Qty/box
AC3930A	Angel [®] Catheter	5
AK9035A	Percutaneous Access Kit	5

MM 02-6000-02 USA Rev. 03.

