



weLLgo™ Chiba Biopsy Needle for Fine Needle Aspiration


The weLLgo™ chiba biopsy needle features a chiba needle tip design to accommodate differences in user preference.

The chiba needle is highly versatile and can be used for fine needle aspiration biopsy and injection of drugs and contrast agents. It also features color coding on the hub of the needle to make identification of accurate gauge size easy.

Features:

- Sharp chiba needle tip designed for minimal tissue damage.
- Adjustable needle stopper for increased penetration depth control.
- Universal needle design for variances in applications to ease procurement and inventory control.

Color-coded hub to make identification of accurate gauge size easy.



Echogenic needle tip design which provides excellent visualization during procedure.

Centimeter markings on cannula provide depth insertion reference.

Ordering Information

weLLgo™ Chiba Biopsy Needle Fine Needle Aspiration

Ref no.	Gauge	Length	Qty/box
061810	18	10 CM	25
061815	18	15 CM	25
061820	18	20 CM	25
061825	18	25 CM	25
062010	20	10 CM	25
062015	20	15 CM	25
062020	20	20 CM	25
062025	20	25 CM	25
062110	21	10 CM	25
062115	21	15 CM	25
062120	21	20 CM	25
062125	21	25 CM	25
062210	22	10 CM	25
062215	22	15 CM	25
062220	22	20 CM	25
062225	22	25 CM	25

Contact Details

www.mermaidmedical.com

customer@mermaidmedical.com

info@mermaidmedical.nl

customer@mermaidmedical.co.uk